

## Technical Bulletin

<b>Product:</b>	2WIN-S – Kaleidos APP
<b>Date:</b>	2023/10/27 (October 27 <sup>th</sup> , 2023)
<b>Bulletin number and title:</b>	34 - New Kaleidos APP release 1.4.0.0
<b>Issue has a relevance for the device safety:</b>	NO
<b>Field Safety Corrective Action:</b>	Not required
<b>Issue has a relevance for the device performance:</b>	YES
<b>Reason:</b>	New Kaleidos APP software version
<b>S/N affected:</b>	NA
<b>Scope of notice:</b>	All Distributors/Other
<b>Service time:</b>	Less than 1 minutes for software update procedure
<b>Service parts:</b>	Google Play Store (App) Installation file: <i>Kaleidos_1.4.0.0_(50)_release.apk</i>
<b>Description:</b>	See next page/s
<b>Required action:</b>	Update recommended
<b>Hardware compatibility:</b>	2WIN: from BOM 29.1; compatible tablets cited in the Operator Manual
<b>Software compatibility:</b>	2WIN: from 5.4.2; tablet: Android 8, 9, 10, 11, 12, 13
<b>For further info please contact:</b>	<p><b>Adaptica Srl</b>  Via San Marco, 9/h – 35129, Padova (PD), Italy  Tel: +39 049 773 968  Fax: +39 049 097 0901  Email: support@adaptica.com</p>

---

## **Summary of changes and bugs fixed in Kaleidos APP SW version 1.4.0.0**

- Fixed minor bugs
- Added Android compatibility: from 8.0 to 13
- Added screening mode