

Technical Bulletin

Product:	2WIN
Date:	2020/07/27 (July 27, 2020)
Bulletin number and title:	27 – New software release 5.5.1
Issue has a relevance for the device safety:	NO
Field Safety Corrective Action:	Not required
Issue has a relevance for the device performance:	YES
Reason:	<ul style="list-style-type: none"> - New 2WIN software version (bug fixing) - New KALEIDOS App version (bug fixing) - New Operator Manual
S/N affected:	All
Scope of notice:	All Distributors/Other
Service time:	Less than 5 minutes for software update procedure
Service parts:	full-system-5.5.1
Description:	See next pages
Required action:	Update recommended
Hardware compatibility:	All
For further info please contact:	<p>Adaptica srl Via San Marco, 9/h – 35129, Padova (PD), Italy Tel: +39 049 773 968 Fax: +39 049 097 0901 Email: support@adaptica.com</p>

Summary of bugs fixed in SW version 5.5.1

- New update for 2WIN software, which fixes “Pupils too small” and “Unreliable data” bugs.
- New update for KALEIDOS application, which fixes the bug in the PDF report creation on Android 9 and Android 10 tablets.
- New Operator Manual.

IMPORTANT NOTE: it is not possible to downgrade 2WIN’s software to a previous version.