

Technical Bulletin

| | |
|--|--|
| Product: | 2WIN-S – Kaleidos APP |
| Date: | 2023/01/09 (January 09 th , 2023) |
| Bulletin number and title: | 33 - New Kaleidos APP release 1.3.13.1 |
| Issue has a relevance for the device safety: | NO |
| Field Safety Corrective Action: | Not required |
| Issue has a relevance for the device performance: | YES |
| Reason: | New Kaleidos APP software version |
| S/N affected: | NA |
| Scope of notice: | All Distributors/Other |
| Service time: | Less than 3 minutes for software update procedure |
| Service parts: | Google Play Store (App) Installation file: <i>2winbt_v45_1.3.13.1-release.apk</i> |
| Description: | See next page/s |
| Required action: | Update recommended |
| Hardware compatibility: | 2WIN: from BOM 29.1; compatible tablets cited in the Operator Manual |
| Software compatibility: | 2WIN: from 5.4.2; tablet: Android 8, 9, 10, 11, 12 |
| For further info please contact: | <p>Adaptica Srl Via San Marco, 9/h – 35129, Padova (PD), Italy Tel: +39 049 773 968 Fax: +39 049 097 0901 Email: support@adaptica.com</p> |

Summary of changes and bugs fixed in Kaleidos APP SW version 1.3.13.1

- Fixed minor bugs
- Added Android compatibility: from 8.0 to 12
- Added possibility to download the application from Google Play Store